



Mass Spectrometry Ionisation Techniques

How do we handle different types of sample ?

Types of ionisation techniques


▶ Volatile samples

- Electron Ionisation
- Chemical Ionisation
- GC (and LC) inlets

▶ Non-volatile samples

- Fast Atom Bombardment
- Thermospray
- Matrix Assisted Laser Desorption Ionisation
- Electrospray Ionisation
- Atmospheric Pressure Chemical Ionisation
- LC (and GC) inlets

Comparison of Ionisation Techniques



ThermoFinnigan


Electron Ionisation

- Widely used technique when coupled to GC
- Suitable for volatile organic compounds
 - eg hydrocarbons, oils, flavours, fragrances
- Not really coupled to LC today
- Also called electron impact

Electron Ionisation


- Produces $M^{+\cdot}$ radical cation giving molecular weight
- Produces abundant fragment ions
- Library searchable spectra
- Energetic process. A heated filament emits electrons which are accelerated by a potential difference of usually 70eV into the sample chamber. Ionisation of the sample occurs by removal of an electron from the molecule thus generating a positively charged ion with one unpaired electron.

Electron Ionisation




ThermoFinnigan

Electron Ionisation



Fragmentation




Chemical Ionisation

- Development from EI
- Same compound classes as EI
- Gives molecular weight
- Softer ionisation technique
- Produces $M+H^+$ ions or $M - H^-$ ions
- Used to produce more abundant molecular ions when the molecule under investigation fragments using EI

Chemical Ionisation

- ▶ Similar ionisation technique to EI except that a reagent gas is introduced into the chamber in excess of the sample
- ▶ Positive CI uses methane, isobutane or ammonia as reagent gases
- ▶ Negative CI uses methane reagent gas in electron capture mode
- ▶ Ionised reagent gas protonate the sample molecules leaving a neutral reagent gas species
- ▶ Not reproducible from lab to lab, hence no CI libraries.


Chemical Ionisation





ThermoFinnigan

Chemical Ionisation


Primary ions



Secondary ions



Proton donation




ThermoFinnigan

Fast Atom Bombardment

- Used for large compounds with low volatility (eg peptides, proteins, carbohydrates)
- Solid or liquid sample is mixed with a non-volatile matrix (eg glycerol, crown ethers, nitrobenzyl alcohol)
- Immobilised matrix is bombarded with a fast beam of Argon or Xenon atoms. Charged sample ions are ejected from the matrix and extracted into the mass analysers
- Gives $M+H^+$ or $M+Na^+$ ions
- Choosing correct matrix is difficult

FAB Source



ThermoFinnigan


Matrix Assisted Laser Desorption Ionisation

- Similar process to FAB
- Sample is dissolved in matrix which absorbs light from a short pulse of laser of a specific wavelength. The sample becomes ionised and extracted towards the mass analysers
- Coupled to Time of Flight MS
- Not coupled to LC
- High mass range achievable
- Calibrants may be external or included in sample
- Reproducibility issues

Thermospray

- First widely used LC/MS interface
- Flow rates 0.5 - 1.5 ml/min
- Good for polar compounds
- LC eluent containing sample and ammonium acetate is pumped through a heated vaporiser. The jet of vapour contains small charged droplets which evaporate under the heat and vacuum expelling charged ions from the surface
- Produces $M+H^+$ or $M - H^-$ ions
- Not commercially available today

Thermospray Process



ThermoFinnigan

Atmospheric Pressure Ionisation


- ▶ Most important and widely used LC / MS technique
- ▶ API two types
 - Electrospray
 - Atmospheric Pressure Chemical Ionisation
- ▶ > 99% new LC/MS use API source
- ▶ Ionisation takes place outside vacuum region

Atmospheric Pressure Ionisation

- ▶ API coupled to LC or CE or Nanospray
- ▶ Handle wide range of flow rates
- ▶ Produce Intense $M+H^+$ ions
- ▶ Very little fragmentation
 - Need MS/MS for structural information
- ▶ Applicable to wide range of compounds
- ▶ Sample must be in solution

Electrospray

- ▶ Electrospray also known as :
 - Ionspray
 - Nanospray
 - Sonic Spray
 - "Pure" Electrospray
 - ESI, ES, IS




ThermoFinnigan

Electrospray


- Softest ionisation technique
- Best for polar non-volatile compounds (proteins, peptides, nucleic acids, Pharmaceuticals, natural products)
- Coupled to LC at a flow range of 2-1000 $\mu\text{L}/\text{min}$, nanospray (10 nL/min - 2 $\mu\text{L}/\text{min}$)
- Ions are ejected from charged vapour droplets to gas phase producing $M+H^+$ or $M - H^-$ ions
- Can produce multiply charged ions allowing determination of high molecular weight proteins
- Not very tolerant of non-volatile salts

ESI Probe



ThermoFinnigan


Electrospray Process



ThermoFinnigan

APCI

- ▶ Atmospheric Pressure Chemical Ionisation, also known as :
 - APCI
 - Heated nebuliser
 - APcl




ThermoFinnigan

APCI


- Used for wide range polarity of compounds
- HPLC eluent (up to 2ml/min flow rate) is vaporised at up to 600 °C
- The Corona discharge needle ionises solvent molecules. A combination of collisions and charge transfer reactions between the solvent and the analyte results in the transfer of a proton to form either $M+H^+$ or $M-H^-$ ions
- Compounds can thermally degrade
- Multiply charged ions rare
- More tolerant to salts

APCI Probe



ThermoFinnigan


APCI Process



Corona discharge needle ionises solvent to generate a chemical ionisation reagent gas plasma

ThermoFinnigan

API Publications



Halket JM and Down S, LC/MS Update, HD Science, Nottingham

ThermoFinnigan

Solvent suitability


■ HPLC buffers

- Reversed phase most often used
- MeOH, ACN, H₂O,
- TFA, formic acid, acetic acid, Ammonium formate, ammonium acetate
- Normal phase can be used

■ Non-volatile buffers


- OSA, aQa self cleaning source, off-axis probe

Orthogonal Sampling Adaptor




ThermoFinnigan

Orthogonal Sampling Adaptor



ThermoFinnigan

OSA fitted - 10mM Phosphate solution




ThermoFinnigan

Problems

How would you analyse this compound ?

Naphthalene




What sample introduction technique could you use ?
Which ionisation technique ?

A: EI, GC/MS

Problems

How would you analyse this compound ?

Phenacetin



What sample introduction technique could you use ?
Which ionisation technique ?

A: API (either APCI or ESI), LC/MS

Problems

How would you analyse myoglobin ?


Myoglobin is a protein with a molecular weight of 16,951.

If the Mass Spectrometer has a mass range of up to 4,000, how can you analyse high molecular weight proteins ?

Multiply charged myoglobin ions from ESI


$$(M_2 - 1.008) / M_1 - M_2 = Z_1$$

$$(Z_1 * M_1) - (Z * 1.008) = M_{wt}$$



ThermoFinnigan

Deconvoluted myoglobin spectrum



ThermoFinnigan